

GSBIO PCR COMPATIBILITY LIST

8-strip Tubes with Attachde Caps												
CAT No.	CP0100	CP0101	CP0110	CP0111	CP0200	CP0201	CP0210	CP0211	CP4000	CP4001	CP4010	CP4011
Applied biosystems/life technologies												
StepOnePlus™ qPCR system	√	√			√	√			√	√		
StepOne™ qPCR system (48孔)	√	√			√	√						
7000 ABI PRISM 7000 SDS(普通扩增)			√	√			√	√			√	√
7300 qPCR system			√	√			√	√			√	√
7500 qPCR system			√	√			√	√			√	√
7500 Fast qPCR system	√	√			√	√			√	√		
7900HT qPCR system			√	√			√	√			√	√
7900HT Fast qPCR system	√	√			√	√			√	√		
9700(普通扩增)												
ViiA™7 fast qPCR 96*0.1ml	√	√			√	√			√	√		
ViiA™7 qPCR system			√	√			√	√			√	√
QuantStudio™ 3/5/6/7/12k Flex Fast	√	√			√	√			√	√		
QuantStudio™ 3/5/6/7/12k Flex			√	√			√	√			√	√
基因测序310/3130/3500/3730												
Bio-Rad Laboratories												
CFX96 Real-Time PCR	√	√			√	√			√	√		
iQ5 Real-Time PCR			√	√			√	√			√	√
iQ Real-Time PCR			√	√			√	√			√	√
MyiQ Real-Time PCR			√	√			√	√			√	√
Option 2	√	√			√	√			√	√		
CFX connect	√	√			√	√			√	√		
CFX 96 DEEP WELL			√	√			√	√			√	√
Roche												
Cobas z 480		√				√				√		
LightCycler® 480 System 384-well		√				√				√		
LightCycler® 480 System		√				√				√		
LightCycler® 96 System		√				√				√		
LightCycler® Nano		√				√				√		
Agilent/Stratagena												
Mx3000P™			√	√			√	√				
Mx3005P™			√	√			√	√				
Mx4000™			√	√			√	√				
ARIA MX G8830A	√	√			√	√			√	√		
Eppendorf												
Mastercycler ep realplex Silver			√	√			√	√			√	√
宏石												
SLAN-48P			√				√				√	
SLAN-96P			√				√				√	
SLAN-96S			√				√				√	
Tianlong Science 天隆												
TL988	√	√			√	√			√	√		
Gentier 48S							√				√	
Gentier 48E							√				√	
Gentier 96E							√				√	
Bioer 博日												
Line-Gene I							√				√	
Line-Gene K 48孔							√				√	
Line-Gene 3310/3320							√				√	
Line-Gene 9600							√				√	
雅睿												
MA6000							√				√	
MA-1600Q扩增仪 (16孔)							√					
MA-3200Q扩增仪 (32孔)							√					
MA-688 (48孔)							√				√	
鲲鹏												
ArchimedX4	√	√			√	√			√	√		

GSBIO PCR COMPATIBILITY LIST

96-well Plates and Assemblies													
CAT No.	CP1000	CP1001	CP1010	CP1011	CP2000	CP2001	CP2010	CP2011	CP2110	CP5000	CP5001	CP5010	CP5011
Applied biosystems/life technologies													
StepOnePlus™ qPCR system	√	√								√	√		
StepOne™ qPCR system (48孔)													
7000 ABI PRISM 7000 SDS(普通扩增)			√	√			√	√	√			√	√
7300 qPCR system			√	√			√	√	√			√	√
7500 qPCR system			√	√			√	√	√			√	√
7500 Fast qPCR system	√	√								√	√		
7900HT qPCR system							√	√	√				
7900HT Fast qPCR system	√	√											
9700(普通扩增)			√				√		√				
ViiA™7 fast qPCR 96*0.1ml	√	√								√	√		
ViiA™7 qPCR system			√	√						√	√		
QuantStudio™ 3/5/6/7/12k Flex Fast	√	√								√	√		
QuantStudio™ 3/5/6/7/12k Flex			√	√								√	√
基因测序310/3130/3500/3730									√				
Bio-Rad Laboratories													
CFX96 Real-Time PCR	√	√			√	√							
iQ5 Real-Time PCR								√					
iQ Real-Time PCR								√					
MyiQ Real-Time PCR								√					
Option 2	√	√			√	√							
CFX connect	√	√			√	√							
CFX 96 DEEP WELL			√	√			√	√	√				
Roche													
Cobas z 480						√							
LightCycler® 480 System 384-well						√							
LightCycler® 480 System						√							
LightCycler® 96 System						√							
LightCycler® Nano						√							
Agilent/Stratagena													
Mx3000P™			√	√									
Mx3005P™			√	√									
Mx4000™			√	√									
ARIA MX G8830A	√	√											
Eppendorf													
Mastercycler ep realplex Silver			√	√					√				
宏石													
SLAN-48P			√				√	√					
SLAN-96P			√				√	√					
SLAN-96S			√				√	√					
Tianlong Science 天隆													
TL988	√	√											
Gentier 48S													
Gentier 48E													
Gentier 96E													
Bioer 博日													
Line-Gene I			√										
Line-Gene K 48孔			√										
Line-Gene 3310/3320			√										
Line-Gene 9600			√										
雅睿													
MA6000			√				√		√				
MA-1600Q扩增仪 (16孔)													
MA-3200Q扩增仪 (32孔)													
MA-688 (48孔)													
鲲鹏													
ArchimedX4	√	√			√	√							